

# Orbitrap Exploris™ Isotope Solutions

Taylor Graham, Brett Davidheiser, Nils Kuhlbusch, Andreas Hilkert, Dieter Juchelka


June 7, 2026

 The world leader in serving science



# Outline

- Magnetic sector portfolio
- Ideal molecule to benchmark orbitrap isotope ratios
- Measurement technique
- Review of established compounds
- Current labs with orbitrap isotope equipment
- Thermo demo lab
- Questions (feel free to ask during presentation)



# The current Thermo Scientific™ IRMS portfolio

We currently offer a scaled portfolio of MS and sample introduction solutions



- **Thermo Scientific™ DELTA Q™ IRMS**
- Sensitivity up to 800 M/I
- Mass range up to  $m/z$  96
- Accommodates up to 10 collectors



- **Thermo Scientific™ 253 Plus™ 10 kV IRMS**
- High-precision isotope analysis
- Long term stability and robustness
- Complete automation for ease-of-use



- **Thermo Scientific™ Orbitrap Exploris™ 120/240/480 MS**
- ESI source
- Variable resolution for multiple isotopologue collection
- Complete automation for ease-of-use

# Oxyanions: refine understanding of element cycling

Atmosphere  $\longleftrightarrow$  Soil  $\longleftrightarrow$  Water




Oxyanions are important chemical constituents of virtually every environment on Earth.

# What are current isotope ratio techniques missing?

## A wish list


- Direct measurement of intact molecules
- Access to liquid, polar samples
- Higher sensitivity
- Less sample preparation, eliminate solid → gas process
- Position specific isotope analysis



Nitrogen **2021**, 2, 196–217. <https://doi.org/10.3390/nitrogen2020013>

# An opportunity to bridge a gap exists


Combine technology developed for Life Science with approaches from stable isotope geochemistry




# Mass spectrometers for isotope ratio analysis

## Classic isotope-ratio MS $N_2, O_2, N_2O$


- **Indirect:** Conversion into **simple gases**
- Combining isotopologues on **few signals**



## Thermo Scientific Orbitrap™ MS



- **Direct:** Intact isotopologues
- Separation of **all species** by HRMS



# Orbitrap Exploris Isotope Solutions


Thermo Scientific™ Orbitrap Exploris™ Isotope Solutions includes:

Thermo Scientific™  
Orbitrap Exploris™  
120/240/480 MS



Dual Syringe Inlet

Data evaluation Package  
for Isotope Ratio MS



Optional Thermo  
Scientific™ Vanquish™  
Neo UHPLC System



# What's inside the Orbitrap Exploris MS

Essential components for isotope ratio analysis




# Important parameters for Orbital analyzer

## Isotope ratio analysis


- Resolution:
  - Orbitrap analyzes a specific number of ions (ion package) per scan
  - Every scan results in a mass spectrum
  - Higher resolution requires longer scan time:

	Resolution at $m/z$ 200	Scans per second
$\text{NO}_3^-$	15,000	22
	30,000	12
$\text{HSO}_4^-$	45,000	10
	60,000	7
Organics	120,000	3
	240,000	1.5



# Isotope ratio methodology for Orbitrap MS


- **Precision: maximizing ion counts**
  - Longer analysis times achieve greater precision.
  - Lower resolution setting requires less time for scanning.
  - Optimize mass range maximize the number of target ions.
- **Accuracy: Sample/Standard comparison**
  - $\text{NO}_3^-$  Analysis of signals (reference/sample) for 6–9 mins



# Dual Syringe Inlet

## Sample introduction technique


- **Dual Syringe Inlet System** utilizing the diverter valve option of the Orbitrap Exploris MS
- Direct infusion of sample (**50  $\mu\text{M}$**  in MeOH) with a flow rate of **4  $\mu\text{l}/\text{min}$**  via a syringe pump; sample reference comparison by switching of a diverter valve.



# Automated In-Flow Injection

## Sample introduction technique


- **In-flow Injection** - Vanquish Neo UHPLC System
- Loop-injection of **20-30  $\mu\text{L}$**  sample (**50  $\mu\text{M}$**  in MeOH) by the Autosampler into a flow of **4  $\mu\text{L}/\text{min}$**  of MeOH.



# Automated In-Flow Injection with Syringe Reference



## Sample introduction technique

- **In-flow Injection** - Vanquish Neo UHPLC System
- Syringe injected reference material to control for ESI instrument drift. Increases orbitrap uptime while loop is washed out.




# Orbitrap for isotopes: workflow for nitrate

Nitrate simplest model (4 atoms with isotope species:  $^{14}\text{N}$ ,  $^{15}\text{N}$ ,  $^{16}\text{O}$ ,  $^{17}\text{O}$ ,  $^{18}\text{O}$ )



# ESI-Orbitrap for isotopes – methodology

## 'M0' experiment




- M0 peak as basepeak

	<i>m/z</i>	Isotopologue	Abundance	Percentage
<b>M0</b>	<b>61.9884</b>	<sup>14</sup> N <sup>16</sup> O <sub>3</sub>	<b>989242</b>	<b>98.9</b>
<b>M+1</b>	<b>62.9854</b>	<sup>15</sup> N <sup>16</sup> O <sub>3</sub>	<b>3637</b>	<b>0.36</b>
	<b>62.9926</b>	<sup>14</sup> N <sup>17</sup> O <sup>16</sup> O <sub>2</sub>	<b>1127</b>	<b>0.11</b>
<b>M+2</b>	63.9896	<sup>15</sup> N <sup>17</sup> O <sup>16</sup> O <sub>2</sub>	4.1	
	<b>63.9926</b>	<sup>14</sup> N <sup>18</sup> O <sup>16</sup> O <sub>2</sub>	<b>5951</b>	<b>0.60</b>
	63.9968	<sup>14</sup> N <sup>17</sup> O <sub>2</sub> <sup>16</sup> O	0.4	
<b>M+3</b>	64.9897	<sup>14</sup> N <sup>18</sup> O <sup>16</sup> O <sub>2</sub>	21.9	
	64.9938	<sup>15</sup> N <sup>17</sup> O <sub>2</sub> <sup>16</sup> O	< 0.1	
	64.9968	<sup>15</sup> N <sup>17</sup> O <sup>18</sup> O <sup>16</sup> O	4.5	
	65.0010	<sup>14</sup> N <sup>17</sup> O <sub>3</sub>	< 0.1	
	<b>M+4</b>	65.9939	<sup>15</sup> N <sup>18</sup> O <sup>16</sup> O <sub>2</sub>	< 0.1
	65.9969	<sup>15</sup> N <sup>18</sup> O <sub>2</sub> <sup>16</sup> O	11.9	
	65.9981	<sup>15</sup> N <sup>17</sup> O <sub>3</sub>	< 0.1	
	66.0011	<sup>15</sup> N <sup>17</sup> O <sub>2</sub> <sup>18</sup> O	< 0.1	

# ESI-Orbitrap for isotopes – methodology

## 'noM0' experiment



- $^{15}\text{N}$  or  $^{18}\text{O}$  peak as base peak

	<i>m/z</i>	Isotopologue	Abundance	Percentage
<del>M0</del>	<del>61.9884</del>	<del><math>^{14}\text{N}^{16}\text{O}_3</math></del>	<del>989242</del>	<del>98.9</del>
<b>M+1</b>	<b>62.9854</b>	<b><math>^{15}\text{N}^{16}\text{O}_3</math></b>	<b>3637</b>	<b>33.81</b>
	<b>62.9926</b>	$^{14}\text{N}^{17}\text{O}^{16}\text{O}_2$	<b>1127</b>	<b>10.48</b>
<b>M+2</b>	63.9896	$^{15}\text{N}^{17}\text{O}^{16}\text{O}_2$	4.1	
	<b>63.9926</b>	<b><math>^{14}\text{N}^{18}\text{O}^{16}\text{O}_2</math></b>	<b>5951</b>	<b>55.32</b>
	63.9968	$^{14}\text{N}^{17}\text{O}_2^{16}\text{O}$	0.4	
<b>M+3</b>	<b>64.9897</b>	<b><math>^{15}\text{N}^{18}\text{O}^{16}\text{O}_2</math></b>	<b>21.9</b>	<b>0.20</b>
	64.9938	$^{15}\text{N}^{17}\text{O}_2^{16}\text{O}$	< 0.1	
	<b>64.9968</b>	<b><math>^{14}\text{N}^{17}\text{O}^{18}\text{O}^{16}\text{O}</math></b>	<b>4.5</b>	<b>0.04</b>
	65.0010	$^{14}\text{N}^{17}\text{O}_3$	< 0.1	
<b>M+4</b>	65.9939	$^{15}\text{N}^{18}\text{O}^{16}\text{O}_2$	< 0.1	
	<b>65.9969</b>	<b><math>^{14}\text{N}^{18}\text{O}_2^{16}\text{O}</math></b>	<b>11.9</b>	<b>0.11</b>
	65.9981	$^{15}\text{N}^{17}\text{O}_3$	< 0.1	
	66.0011	$^{14}\text{N}^{17}\text{O}_2^{18}\text{O}$	< 0.1	


# Applications

Published or application notes available for:

- Nitrate
- Sulfate
- Phosphate
- MSA
- Acetate
- Amino Acids – glutamic acid and alanine
- Caffeine
- Vanillin
- Methyl Phosphonic acid (chemical weapon forensics)
- DDT
- PFAS
- SMX
- Glucose
- FAME
- Neonicotinoids
- Perchlorate

# Questions

1. What is the ionization method for the Orbitrap?
2. Where can we fragment our molecule?



# Dual Inlet measurements of nitrate


## 'M0' experiment:

- Three nitrate reference materials available from USGS:[1]

	$\delta^{15}\text{N}_{\text{AIR}}$	$\delta^{18}\text{O}_{\text{VSMOW}}$	$^{17}\text{O}_{\text{VSMOW}}$
USGS32	+ 180 ‰	+ 25.7 ‰	-
USGS34	- 1.8 ‰	- 27.9 ‰	- 14.8 ‰
USGS35	+2.7 ‰	+ 57.5 ‰	+ 51.5 ‰

## Quality control: $^{15}\text{N}/^{14}\text{N}$ isotope ratio analysis

Block	Description	Ratio ( $^{15}\text{N}/^{14}\text{N}$ )	$\delta^{15}\text{N}_{\text{USGS32}/\text{Air}}$ [‰]	$\delta^{15}\text{N}_{\text{USGS32}/\text{Air}}$ [‰]	Std. Dev.
1	Reference (USGS35)	0.00430	179.5	179.6	0.4
2	Sample (USGS32)	0.00506			
3	Reference (USGS35)	0.00430			
4	Sample (USGS32)	0.00506	179.3		
5	Reference (USGS35)	0.00430	180.1		
6	Sample (USGS32)	0.00506			
7	Reference (USGS35)	0.00430			



[1] Böhlke, J.K. ; et al. Rapid Commun. Mass Spectrom., 2003, 17, p. 1835–1846.


# Referencing scheme for nitrate

## Analysis of international standards for nitrate

- Isotope ratio data of 'M0' and 'noM0' experiments using N11 as a working standard:[1]

	$\delta^{15}\text{N}_{\text{AIR}}$	$\delta^{18}\text{O}_{\text{VSMOW}}$	$^{17}\text{O}_{\text{VSMOW}}$
USGS32	+ 180 ‰	+ 25.7 ‰	-
USGS34	- 1.8 ‰	- 27.9 ‰	- 14.8 ‰
USGS35	+2.7 ‰	+ 57.5 ‰	+ 51.5 ‰


Isotope ratio data of six of nitrate's most abundant isotopologues using the In-Flow Injection setup.



[1] Böhlke, J.K. ; et al. Rapid Commun. Mass Spectrom., 2003, 17, p. 1835–1846.

# Nitrate analysis of environmental samples

## 'M0' experiment





## Dual Syringe Inlet methodology

- Sample prep: "Dilute and Shoot"
  - Dilution to **1 $\mu\text{M}$  (62 pg  $\text{NO}_3^-/\mu\text{L}$ )** for high salinity samples
- Results  $\delta^{15}\text{N}$  and  $\delta^{18}\text{O}$ 
  - Precision: < 1‰
  - Difference Orbitrap – Expected: < 1‰
- Takeaway
  - ESI tolerates up to **100-fold  $\text{Cl}^-$**  load
  - Ground waters diluted up to **1/500**

# Other oxyanions – sulfate and phosphate

Sulfate ionized as  $\text{HSO}_4^-$

Phosphate ionized as  $\text{H}_2\text{PO}_4^-$



## Stable Isotope Analysis of Intact Oxyanions Using Electrospray Quadrupole-Orbitrap Mass Spectrometry



Cajetan Neubauer,\* Antoine Crémère, Xingchen T. Wang, Nivedita Thiagarajan, Alex L. Sessions, Jess F. Adkins, Nathan F. Dalleska, Alexandra V. Turchyn, Josephine A. Clegg, Annie Moradian, Michael J. Sweredoski, Spiros D. Garbis, and John M. Eiler

Cite This: *Anal. Chem.* 2020, 92, 3077–3085

Read Online

# Isotope mapping of sulfates

	<i>m/z</i>	Isotopologue	Abundance
<b>M0</b>	<b>96.9596</b>	$1\text{H}^{32}\text{S}^{16}\text{O}_4$	9408592
<b>M+1</b>	<b>97.9590</b>	$1\text{H}^{33}\text{S}^{16}\text{O}_4$	7427
	<b>97.9638</b>	$1\text{H}^{32}\text{S}^{17}\text{O}^{16}\text{O}_3$	1433
<b>M+2</b>	<b>98.9554</b>	$1\text{H}^{34}\text{S}^{16}\text{O}_4$	42084
	<b>98.9638</b>	$1\text{H}^{32}\text{S}^{18}\text{O}^{16}\text{O}_3$	7632
<b>M+3</b>	<b>99.9596</b>	$1\text{H}^{34}\text{S}^{17}\text{O}^{16}\text{O}_3$	63
	<b>99.9632</b>	$1\text{H}^{33}\text{S}^{18}\text{O}^{16}\text{O}_3$	59
	<b>99.9680</b>	$1\text{H}^{32}\text{S}^{17}\text{O}^{18}\text{O}^{16}\text{O}_2$	9
<b>M+4</b>	<b>100.9546</b>	$1\text{H}^{36}\text{S}^{16}\text{O}_4$	145
	<b>100.9596</b>	$1\text{H}^{34}\text{S}^{18}\text{O}^{16}\text{O}_3$	333
	<b>100.9680</b>	$1\text{H}^{32}\text{S}^{18}\text{O}_2^{16}\text{O}_2$	44



# Isotopes of Phosphate





Figure 2. Proportional variations of three ions, including PO<sub>2</sub><sup>-</sup> (P<sup>16</sup>O<sub>2</sub><sup>-</sup>  $m/z$  = 63, P<sup>17</sup>O<sup>16</sup>O<sup>-</sup>  $m/z$  = 64, P<sup>18</sup>O<sup>16</sup>O<sup>-</sup>  $m/z$  = 65), PO<sub>3</sub><sup>-</sup> (P<sup>16</sup>O<sub>3</sub><sup>-</sup>  $m/z$  = 79, P<sup>17</sup>O<sup>16</sup>O<sub>2</sub><sup>-</sup>  $m/z$  = 80, P<sup>18</sup>O<sup>16</sup>O<sub>2</sub><sup>-</sup>  $m/z$  = 81), and H<sub>2</sub>PO<sub>4</sub><sup>-</sup> (<sup>1</sup>H<sub>2</sub>P<sup>16</sup>O<sub>4</sub><sup>-</sup>  $m/z$  = 97, <sup>1</sup>H<sub>2</sub>P<sup>17</sup>O<sup>16</sup>O<sub>3</sub><sup>-</sup>  $m/z$  = 98, <sup>1</sup>H<sub>2</sub>P<sup>18</sup>O<sup>16</sup>O<sub>3</sub><sup>-</sup>  $m/z$  = 99) observed at different collision energies in both PO<sub>3</sub><sup>-</sup> out of source fragmentation (A) and PO<sub>3</sub><sup>-</sup> out of HCD (B).

analytical  
chemistry

pubs.acs.org/ac

Article




## Oxygen Isotope Analysis of Nanomole Phosphate Using PO<sub>3</sub><sup>-</sup> Fragment in ESI-Orbitrap-MS

Zhenfei Wang, Shohei Hattori, Yongbo Peng,\* Longchen Zhu, Zhao Wei, and Huiming Bao

Cite This: <https://doi.org/10.1021/acs.analchem.3c03070>

Read Online


# First results of methanesulfonic acid (MSA)



# Fragmentation of methanesulfonic acid (MSA)


240k Resolution

~ 1 scan/s



45k Resolution

~ 10 scan/s



# MSA using $\text{SO}_3^-$ Fragment Paper

pubs.acs.org/ac

Article

## Triple Oxygen Isotope Analysis of Methanesulfonate Using the $\text{SO}_3^-$ Fragment in ESI-Orbitrap-MS

Yihang Hong, Longchen Zhu, Daniel R. Crocker, Tengyu Liu, Zhenfei Wang, Zhao Wei, Chen Yu, Yu Wei, Hao Yan, David T. Johnston, Issaku E. Kohl, Yongbo Peng, Andreas Hilkert, Cajetan Neubauer, and Shohei Hattori\*



Cite This: *Anal. Chem.* 2025, 97, 14339–14348



Read Online

- HCD Collision Energy ~60% (compared to 50% for  $\text{PO}_3^-$ )
- Offline purification and separation to remove  $\text{SO}_3^-$  contribution to MSA  $\text{SO}_3^-$
- Murchison meteorite contains distinctly extraterrestrial isotope signatures for  $\delta^{13}\text{C}$ ,  $\delta\text{D}$  and  $\Delta^{33}\text{S}$





Figure 2. Proportion variations of three ions of  $\text{SO}_2^-$  [ $^{32}\text{S}^{16}\text{O}_2^-$  ( $m/z = 63.96245$ )],  $\text{SO}_3^-$  [ $^{32}\text{S}^{16}\text{O}_3^-$  ( $m/z = 79.95736$ )], and  $\text{CH}_3\text{SO}_3^-$  [ $^{12}\text{C}^1\text{H}_3^{32}\text{S}^{16}\text{O}_3^-$  ( $m/z = 94.98084$ )] under the different (A) in-source fragmentation energy and (B) HCD collision energy conditions.


# Acetate Isotopes on the Orbitrap



Analytical Chemistry

pubs.acs.org/ac

Article



Article

## Simultaneous, High-Precision Measurements of $\delta^2\text{H}$ and $\delta^{13}\text{C}$ in Nanomole Quantities of Acetate Using Electrospray Ionization-Quadrupole-Orbitrap Mass Spectrometry

Elliott P. Mueller,\* Alex L. Sessions, Peter E. Sauer, Gabriella M. Weiss, and John M. Eiler



Cite This: *Anal. Chem.* 2022, 94, 1092–1100



Read Online

# Acetate Isotopes on the Orbitrap





Figure 3: Five mixtures of organic acids containing different proportions of acetate, propionate, and butyrate with known  $\delta^{13}\text{C}$  and  $\delta^2\text{H}$  values were measured against a pure acetate standard. All mixtures had an absolute acetate concentration of 50  $\mu\text{M}$ , yet acetate's SIC changed dramatically.



[pubs.acs.org/ac](https://pubs.acs.org/ac)

Article

## Simultaneous, High-Precision Measurements of $\delta^2\text{H}$ and $\delta^{13}\text{C}$ in Nanomole Quantities of Acetate Using Electrospray Ionization-Quadrupole-Orbitrap Mass Spectrometry

Elliott P. Mueller,\* Alex L. Sessions, Peter E. Sauer, Gabriella M. Weiss, and John M. Eiler



Cite This: *Anal. Chem.* 2022, 94, 1092–1100



Read Online

# Glutamic Acid Isotopes





Figure 2. MS2 mass spectrum with 70% NCE acquired with glutamic acid.






Figure 3. Isotopic pattern of fragment at m/z 84 at 90,000 resolution.

# Glutamic Acid Isotopes



Thermo Fisher Scientific TN003099

**Table 1**

Details about different methods of position-specific isotope analysis (PSIA). Strengths and weaknesses of each type of PSIA are presented. It is important to note that individual analytical goals and available resources are the overall determining factor for which method is best for a given project.


Method	Measured species	Sensitivity (nmol)	Strengths	Weaknesses	Select References
Chemical or enzymatic degradation + IRMS	Conversion of analyte to CO <sub>2</sub>	100s	Uses established IRMS methods, ability to constrain all positions	Many steps, requires derivatization, large sample size, long analysis times	[2–4]
Pyrolysis + IRMS	Conversion of analyte to CO <sub>2</sub>	100s	Uses established IRMS methods, ability to constrain all positions	Many steps, requires derivatization, large sample size, long analysis times	[5–7]; Dias et al. (2022)
NMR	Intact molecule	50,000–300,000	Ability to constrain all positions, non-destructive, derivatization not required	large sample size (100s of mmol)	[9,10]; Rasmussen and Hoffman (2020)
Orbitrap-MS	Molecular ion, fragment ions	1–10	Measures intact molecule, small sample size, short analysis time, derivatization not required	Currently no on-line separation for ESI-Orbitrap, need for internationally recognized standards, cannot access all positions without chemical preparation	[13,14]; this study



Contents lists available at ScienceDirect

International Journal of Mass Spectrometry

journal homepage: [www.elsevier.com/locate/ijms](http://www.elsevier.com/locate/ijms)



Analysis of intramolecular carbon isotope distributions in alanine by electrospray ionization Orbitrap mass spectrometry

Gabriella M. Weiss<sup>a,b,\*</sup>, Alex L. Sessions<sup>a</sup>, Maxime Julien<sup>c,d</sup>, Timothy Csernica<sup>e</sup>, Keita Yamada<sup>d</sup>, Alexis Gilbert<sup>d,f</sup>, Katherine H. Freeman<sup>b</sup>, John M. Eiler<sup>a</sup>

<sup>a</sup> Division of Geological and Planetary Sciences, California Institute of Technology, Pasadena, CA, 91125, USA

<sup>b</sup> Department of Geosciences, The Pennsylvania State University, University Park, PA, 16803, USA




<sup>c</sup> GFZ German Research Centre for Geosciences - Helmholtz Centre Potsdam, Organic Geochemistry, Potsdam, Germany

<sup>d</sup> Department of Chemical Science and Engineering, Tokyo Institute of Technology, Yokohama, 226-8503, Kanagawa, Japan

<sup>e</sup> Division of Chemistry and Chemical Engineering, California Institute of Technology, Pasadena, CA, 91125, USA

<sup>f</sup> Earth-Life Science Institute, Tokyo Institute of Technology, Meguro, 152-8550, Tokyo, Japan

# Alanine Isotopes




ELSEVIER

Contents lists available at ScienceDirect

International Journal of Mass Spectrometry

journal homepage: [www.elsevier.com/locate/ijms](http://www.elsevier.com/locate/ijms)



Analysis of intramolecular carbon isotope distributions in alanine by electrospray ionization Orbitrap mass spectrometry

Gabriella M. Weiss<sup>a,b,\*</sup>, Alex L. Sessions<sup>a</sup>, Maxime Julien<sup>c,d</sup>, Timothy Csernica<sup>e</sup>, Keita Yamada<sup>d</sup>, Alexis Gilbert<sup>d,f</sup>, Katherine H. Freeman<sup>b</sup>, John M. Eiler<sup>a</sup>

<sup>a</sup> Division of Geological and Planetary Sciences, California Institute of Technology, Pasadena, CA, 91125, USA

<sup>b</sup> Department of Geosciences, The Pennsylvania State University, University Park, PA, 16803, USA



<sup>c</sup> GFZ German Research Centre for Geosciences - Helmholtz Centre Potsdam, Organic Geochemistry, Potsdam, Germany

<sup>d</sup> Department of Chemical Science and Engineering, Tokyo Institute of Technology, Yokohama, 226-8503, Kanagawa, Japan

<sup>e</sup> Division of Chemistry and Chemical Engineering, California Institute of Technology, Pasadena, CA, 91125, USA

<sup>f</sup> Earth-Life Science Institute, Tokyo Institute of Technology, Meguro, 152-8550, Tokyo, Japan

# Accuracy - Caffeine reference materials




	$\delta^{13}\text{C}_{\text{VPDB, reported}}$	$\delta^{13}\text{C}_{\text{VPDB, measured}}$
IAEA600*	-27.8	-27.8
USGS61	-35.1	-34.7 ± 0.7
USGS62	-14.8	-14.1 ± 0.6
USGS63	-1.2	-1.3 ± 0.6
Lab_Std.	-33.2	-31.4 ± 0.6

\*Used as calibration point



- Work in press by Nils Kuhlbusch

# Position specific labels in reference materials




• Work in press by Nils Kuhlbusch

# Position specific isotope analysis (PSIA)




• Work in press by Nils Kuhlbusch


# PSIA of caffeine reference materials



USGS61




USGS63




• Work in press by Nils Kuhlbusch



# Vanillin data



2-point  
calibration



# Methyl Phosphonic Acid



ELSEVIER

Talanta

journal homepage: [www.elsevier.com/locate/talanta](http://www.elsevier.com/locate/talanta)

Simultaneous observation of  $^2H$  and  $^{13}C$  enrichment of methyl phosphonic acid via Orbitrap-IRMS with applications to nerve agent forensics

Timothy Csernica<sup>a,\*</sup>, James J. Moran<sup>b,c</sup>, Carlos G. Fraga<sup>d</sup>, John M. Eiler<sup>a</sup>


<sup>a</sup> Division of Geological and Planetary Sciences, California Institute of Technology, Pasadena, CA, 91125, USA

<sup>b</sup> Department of Integrative Biology and Department of Plant, Soil, and Microbial Sciences, Michigan State University, East Lansing, MI, USA

<sup>c</sup> Pacific Northwest National Laboratory, Richland, WA, USA

<sup>d</sup> Air Force Research Laboratory, 10 E. Saturn Blvd, Edwards Air Force Base, CA, 93524, USA

# Methyl Phosphonic Acid



Talanta

journal homepage: [www.elsevier.com/locate/talanta](http://www.elsevier.com/locate/talanta)

Simultaneous observation of  $^2\text{H}$  and  $^{13}\text{C}$  enrichment of methyl phosphonic acid via Orbitrap-IRMS with applications to nerve agent forensics






Timothy Csernica<sup>a,\*</sup>, James J. Moran<sup>b,c</sup>, Carlos G. Fraga<sup>d</sup>, John M. Eiler<sup>a</sup>

<sup>a</sup> Division of Geological and Planetary Sciences, California Institute of Technology, Pasadena, CA, 91125, USA

<sup>b</sup> Department of Integrative Biology and Department of Plant, Soil, and Microbial Sciences, Michigan State University, East Lansing, MI, USA

<sup>c</sup> Pacific Northwest National Laboratory, Richland, WA, USA

<sup>d</sup> Air Force Research Laboratory, 10 E. Saturn Blvd, Edwards Air Force Base, CA, 93524, USA



Contents lists available at [ScienceDirect](https://www.sciencedirect.com)

Marine Pollution Bulletin

journal homepage: [www.elsevier.com/locate/marpolbul](http://www.elsevier.com/locate/marpolbul)




## Multi-isotope characterization of 4,4'-DDT isolated from sediment samples via Orbitrap-IRMS

Timothy A. Csernica<sup>a,\*</sup>, John M. Eiler<sup>a</sup>, Alex L. Sessions<sup>a</sup>

<sup>a</sup> Division of Geological and Planetary Sciences, California Institute of Technology, Pasadena, CA, 91125, United States of America

# Carbon stable isotope analysis of PFAS sources

- Carbon isotope analysis of PFAS, PFOA and FTOH were acquired from multiple vendors and where possible, from multiple lots.
- There was variation in the carbon isotope values among different compounds, manufacturers, and batches.
  - Ex. PFOA Syn  $\delta^{13}\text{C}$  lots differed ( $-28.9 \pm 0.3\text{‰}$ , and,  $-46.1 \pm 0.7\text{‰}$ ).
- Variability among mean carbon isotope values makes EA-IRMS a potentially powerful tracking tool for PFAS contaminants.



- Dombrowski et al (2025). *Environmental Science & Technology Letters* <https://doi.org/10.1021/acs.estlett.5c00021>



# PFOA Carbon Isotope “Fingerprinting”

RESEARCH ARTICLE **OPEN ACCESS**

## Stable Carbon Isotope Analysis of Perfluorooctanoic Acid (PFOA) by Microflow-High Pressure Liquid Chromatography-Orbitrap Mass Spectrometry

Paul K. Wojtal<sup>1,2</sup>  | Brett Davidheiser-Kroll<sup>3</sup>  | Chad S. Lane<sup>1,2</sup>  | Ralph N. Mead<sup>1,2</sup> 

<sup>1</sup>Center for Marine Science, University of North Carolina Wilmington, Wilmington, North Carolina, USA | <sup>2</sup>Department of Earth and Ocean Sciences, University of North Carolina Wilmington, Wilmington, North Carolina, USA | <sup>3</sup>Thermo Fisher Scientific, Bremen, Germany




## Stable Carbon Isotope Analysis of PFOA and PFOS Using Orbitrap Mass Spectrometry

Hiroto Kawashima,\* Tomoha Iezumi, Momoka Suto, and Sachi Taniyasu

Cite This: *Environ. Sci. Technol. Lett.* 2026, 13, 574–579

Read Online





# Peak broadening with heart-cut setup

- Caffeine reference in MeOH + 0.1 % FA
- Sample injected onto RPLC column
- MeOH + 0.1 % FA flushing the heart-cut loop


Can be:

- Syringe
- Additional pump
- Autosampler





# Peak broadening with heart-cut setup

- Broadening of peaks via heart-cut to maximize the number of ions counted



# Caffeine in beverages and food



# HPLC Isotope Orbitrap Coupling

## Coupling Liquid Chromatography to Orbitrap Isotope Ratio Mass Spectrometry: Overcoming Isotope Effects of Chromatography and Amount-Dependency by Peak Homogenization

Aoife Canavan, Leonhard Precht, Habib Al-Ghoul, Nils Kuhlbusch, Andrea M. Erhardt, and Martin Elsner\*

Cite This: *Anal. Chem.* 2026, 98, 590–600

[Read Online](#)

- Sulfamethoxazole as target compound
- 4 nmol of sample
- Reproducible compared to GC-IRMS within 1.5 and 0.9 of  $\delta^{13}\text{C}$  and  $\delta^{34}\text{S}$ , respectively.




Figure 1. Schematic overview of the coupled HPLC-ESI-Orbitrap-MS system for stable isotope analysis.

# HPLC Isotope Orbitrap Coupling

## Coupling Liquid Chromatography to Orbitrap Isotope Ratio Mass Spectrometry: Overcoming Isotope Effects of Chromatography and Amount-Dependency by Peak Homogenization

Aoife Canavan, Leonhard Precht, Habib Al-Ghoul, Nils Kuhlbusch, Andrea M. Erhardt, and Martin Elsner\*

Cite This: *Anal. Chem.* 2026, 98, 590–600

Read Online

- Sulfamethoxazole as target compound
- 4 nmol of sample
- Reproducible compared to GC-IRMS within 1.5 and 0.9 of  $\delta^{13}\text{C}$  and  $\delta^{34}\text{S}$ , respectively.

Scheme 1. Target Fragments of Sulfamethoxazole after Positive or Negative Electrospray Ionization Using Source Fragmentation or HCD Fragmentation in the IRM





Table 1. Instrument Settings Used for Each Fragment of SMX


	F64	F92	F99
polarity ( $\pm$ )	negative	negative	positive
ion transfer tube temperature ( $^{\circ}\text{C}$ )	320	320	320
nominal Orbitrap resolution	30,000	90,000	90,000
AQT isolation window ( $m/z$ )	249.5–254.5	88.5–95.5	94–104
absolute AGC target (#)	$1.5 \times 10^5$	$1.5 \times 10^5$	$1.5 \times 10^5$
microscans (#)	10	10	10
RF lens (%)	100	100	100
maximum ion injection time (ms)	1000	1000	1000
source fragmentation (yes/no)	no	yes	yes
source fragmentation energy (V)	-	45	40
HCD collision (yes/no)	yes <sup>a</sup>	no	no
absolute HCD collision energy (V)	40	-	-

<sup>a</sup>  $m/z$  252 was selected as a precursor ion.

# PSIA of cellulose and glucose by Orbitrap Mass Spectrometry




Thermo Q Exactive Orbitrap:  
**Fragmentation + High-resolution MS**




Hannah  
PhD Student



# Independent observation of $\delta^{13}\text{C}$ of 5 atomic positions in glucose



# Position-specific $\delta^{13}\text{C}$ analysis of C3 and C4 glucose




C4: Corn



C3: Wheat

# Position-specific $\delta^{13}\text{C}$ analysis of C3 and C4 glucose



- 0.1 to 0.5‰ precision on fragment measurements of natural glucose
- Average 1.4‰ precision on calculated position-specific values
- Accurate to externally verified values within precision
- Consumes < 50  $\mu\text{g}$  glucose!

## Position-Specific Carbon Isotope Analysis of Glucose at Natural Isotope Abundance by Electrospray-Ionization Orbitrap Mass Spectrometry

Hannah Dion-Kirschner,\* Celia Kong-Johnson, Kimberly R. Sharp, Nathan Dalleska, John M. Eiler, and Alex L. Sessions




Cite This: *Anal. Chem.* 2026, 98, 915–926



Read Online

# Unique intramolecular $^{13}\text{C}$ patterns in gymnosperm and angiosperm cellulose



*Larix decidua*  
(conifer, dry-adapted)




*Fagus sylvatica*  
(angiosperm, humid-adapted)

## Advancing Stable Isotope Analysis with Orbitrap-MS for Fatty Acid Methyl Esters and Complex Lipid Matrices

Gabriel F. dos Santos,\* Giovanni B. Bevilaqua, Alexis Gilbert, Hugo G. Machado, Maxime Julien, Gesiane S. Lima, Nerilson M. Lima, Júlio C. O. Ribeiro, Alexandre A. Ferreira, Ygor S. Rocha, and Boniek Gontijo\*

Cite This: <https://doi.org/10.1021/jasms.5c00092>

Read Online



Stearic acid





Figure 2. Comparison of  $\delta^{13}\text{C}$  values for stearic acid methyl ester standards obtained using Dual Inlet and HPLC ESI-Orbitrap MS.  $\delta^{13}\text{C}$  values for  $[\text{M} + \text{H}]^+$  and  $[\text{M} + \text{Na}]^+$  ions, showing deviations from the GC-C-IRMS measurements. Error bars represent standard deviations across replicates. Additionally, a summary of MAEs and average deviations for each ionization pathway are shown.

# Organic Acid Review Paper

analytical  
chemistry

Open Access

This article is licensed under [CC-BY 4.0](https://creativecommons.org/licenses/by/4.0/) 

[pubs.acs.org/ac](https://pubs.acs.org/ac)

Article

## Lessons Learned in Orbitrap MS-Based Isotope Ratio Analysis of Organic Acid Mixtures


Hugo G. Machado,\* Elliott P. Mueller, Júlio C. O. Ribeiro, Giovanni B. Bevilaqua, Gabriel F. dos Santos, Alexandre A. Ferreira, Ygor S. Rocha, Surjyendu Bhattacharjee, John M. Eiler, and Boniek Gontijo\*



Cite This: *Anal. Chem.* 2026, 98, 9764–9775



Read Online



# Neonicotinoids Fingerprinting

## Toward Molecular Forensics of Agrochemicals with Orbitrap IRMS: Isotopic Fingerprints of Imidacloprid Sources and Elucidating Hydrolysis

Felix Niemann, Nils J. Kuhlbusch, Annika Gruhlke, Sarah P. Rockel, Robert G. H. Marks, Klaus Kerpen, Milen Nachev, Maik A. Jochmann,\* Andreas Hilkert, and Torsten C. Schmidt




Cite This: <https://doi.org/10.1021/jasms.5c00444>



Read Online

- Faster than LC-IRMS
- Gives multiple isotope ratios from a single analysis
- LOD is much lower than classical IRMS
- Might be able to build databases for complex molecules to fingerprint these compounds.



## Sub-10 Nanomole Perchlorate $\delta^{37}\text{Cl}$ , $\delta^{18}\text{O}$ , and $\Delta^{17}\text{O}$ Measurements by ESI-Orbitrap-MS

Longchen Zhu, Yihang Hong, Shohei Hattori, Zhenfei Wang, Zhao Wei, Yongbo Peng, Nils Johannes Kuhlbusch, and Huiming Bao\*

Cite This: *Anal. Chem.* 2026, 98, 6034–6044


Read Online

- 4 orders of magnitude less sample than conventional IRMS!


Analytical Chemistry

pubs.acs.org/ac

Article




# OEIS Systems in the US



# Demo Lab at CU Boulder

- Orbitrap 240 with Vanquish Neo UHPLC
- Thermo Fisher: Brett Davidheiser
- PIs: Sebastian Kopf and Cajetan Neubauer
- Post Doc: Kristyna Kantnerova (now Assistant Prof at UCT Prague)



# Thank you to our collaborators



Cajetan Neubauer and Kristýna Kantnerová for close collaboration and help

Joel Savarino, John K. Böhlke and Shohei Hattori for providing sample and reference materials

Thermo Fisher Scientific :

Kyle Fort

Konstantin Aizikov

Alexander Makarov

Uwe Rickens

Tabiwang Arrey

Colin Wirth

Christian Klaas

Charles Cartwright

John K. Böhlke  
Stan Mroczkowski



John Eiler  
Nathan Dalleska



Sebastian Kopf



Xingchen Wang




Ralph Mead  
Chad Lane  
Paul Wojtal



# 2<sup>nd</sup> Annual User Meeting July 11<sup>th</sup> and 12<sup>th</sup> at Goldschmidt



My Goldschmidt



## Second Annual International Orbitrap-IRMS Users Meeting

- Home
  - Program by Day
  - Program Overview
  - Exhibitors
  - Program by Theme
  - Program Highlights
  - Early Career
  - Registration
  - Presenter Information
  - Conference Format
  - General Information
  - My Goldschmidt
  - Sign out
  - Technical Support
- Click to add an item to 'My Schedule'.
- Click to add/remove an item to 'My Favorites'.



Saturday, 11 July 2026

08:30 - 17:00

Sunday, 12 July 2026

08:30 - 17:00

**Theme:** Science Workshops


### Description

The use of Orbitrap-based mass spectrometry for stable isotope analysis ('Orbitrap-IRMS') has grown significantly over the past ~10 years, and we expect that it will only continue to grow. Following upon the success of our inaugural meeting at Goldschmidt 2025, our second meeting at Goldschmidt 2026 aims to continue the development of a more formal Community of Practice among Orbitrap-IRMS users (and those interested in adopting the technology). We thus seek participation from individuals at all career stages and particularly encourage early- and mid-career scientists to attend.

This two-day workshop will consist of a brief introduction to Orbitrap-based mass spectrometry for isotope analysis, followed by oral presentations and posters interleaved with moderated discussion sessions. These sessions will focus primarily on scientific results but will also include topics like new analytical methods, interlaboratory comparisons, data processing, etc. The workshop will culminate with a discussion led by the organizers at the end of the second day. In this 'synthesis discussion', we will revisit and summarize the preceding two days' events, solicit feedback for improvement and suggestions for next year's follow-on meeting, identify community 'needs' and 'wants', and select 1–2 new co-convenors who will work with one of us to organize and plan the following meeting at Goldschmidt 2027. We anticipate that this workshop will lead to new collaborations among not only the current Orbitrap-IRMS users but also those interested in entering the field. Ultimately, we aim to initiate a self-perpetuating, community-driven meeting that recurs annually for the first few years before transitioning to a bi-annual cadence.

Oral presentations will be by invitation only, but those who wish to present original research are encouraged to bring a poster. We welcome contributions that present scientific results, scientific works-in-progress, novel methodological approaches and/or instrument adaptations, and "lessons learned" that touch on these or any other related Orbitrap-

# Thank you



## ESI-Orbitrap-IRMS Analysis of Free Metabolites for Early Breast Cancer Detection and Insight into Metabolic Mechanisms

POSTER

Presented by Simon Andren

 Wednesday, 9 July 2025
 10:30 - 12:30
 Exhibit and Poster Hall (2nd and 3rd Floors, Prague Congress Centre)

**Subsession:** 06i-P1 - Technological innovations and methodological advances in isotope geochemistry and isotope ratio mass spectrometry

**Session:** 06i - Technological innovations and methodological advances in isotope geochemistry and isotope ratio mass spectrometry

**Theme:** Theme 06: Frontiers in Analytical and Computational Techniques


### Abstract

Cancer cells rewire their metabolism as they proliferate, altering pathway activities and metabolite quantities. Isotopic analysis of metabolites in blood plasma could thus work as a non-invasive detection method and, in tissue, provide insight into the cancer metabolic network. Studies have shown cancerous cells to be  $^{15}\text{N}$ -depleted and  $^{13}\text{C}$ -enriched relative to adjacent non-cancerous cells[1]; however, primary isotopic anomalies occur at atomic sites undergoing bond-breaking and are not fully expressed in molecular average or bulk isotope analysis. With the advent of Orbitrap-IRMS, site-specific isotope analysis (SSIA) is now possible, which targets the atomic sites in the individual metabolites that experience the isotopic effect from the metabolically dysregulated reaction.

In this study, we focus on two metabolic hallmarks of breast cancer metabolism: the Warburg effect and glutamine addiction. We aim to determine whether SSIA of  $^{13}\text{C}$  in free metabolites of Lactate and Alanine, and  $^{15}\text{N}$  in free metabolites of Glutamine, and Arginine can serve as reliable biomarkers for breast cancer and whether specific isotopic patterns are unique for certain subtypes of breast cancer. Our approach involves three distinct stages: sample purification using reverse-phase HPLC and IC fraction collection, ESI-orbitrap-IRMS measurements, and data processing following the mathematical framework for isotomics. Preliminary quantification from the HPLC of these metabolites in the tissue samples reveals significant differences in cancerous versus adjacent healthy tissue, with cancerous tissues showing multi-fold higher concentrations, along with inter-cancer subtype variations. Using ESI-Orbitrap-IRMS, we achieve sub-permille accuracy on the molecular average for glutamine and approximately 1 permille accuracy for lactate on pure standards, verified by EA measurements. Future work will focus on sample purification and verifying the accuracy of SSIA before measuring biological samples. This study continues the work of SSIA on small quantities of polar organic molecules present in complex matrices. Reliable measurements of the intra-molecular isotopic fingerprints of these metabolites can provide valuable insight into the metabolic mechanisms that sustain cancer cell survival and can be used as a novel breast cancer diagnostic tool.


[1] Tea, I. *et al.*  $^{13}\text{C}$  and  $^{15}\text{N}$  natural isotope abundance reflects breast cancer cell metabolism. *Sci. Rep.* 6, 34251 (2016)

# PFAS




PFOA

Chemical Formula:  $C_8HF_{15}O_2$   
Exact Mass: 413.97370




PFOS

1763-23-1  
Chemical Formula:  $C_8HF_{17}O_3S$   
Exact Mass: 499.93749




PFBA

Chemical Formula:  $C_4HF_7O_2$   
Exact Mass: 213.98648




TFA

Chemical Formula:  $C_2HF_3O_2$   
Exact Mass: 113.99286




PFMOAA 674-13-5

Chemical Formula:  $C_3HF_5O_3$   
Exact Mass: 179.98458




HFPO-DA  
13252-13-6

Chemical Formula:  $C_3F_7O^*$   
Exact Mass: 184.98374



Chemical Formula:  $C_3HF_4O_2^*$   
Exact Mass: 144.99127

HFPO-DA in-source fragmentation very common




PFBS

375-73-5  
Chemical Formula:  $C_4HF_9O_3S$   
Exact Mass: 299.95027

# Outlook: Position Specific Isotope Analysis (PSIA)


## Vanillin

- Fragmentation pattern<sup>[6]</sup> defines PSIA targets
- Dual Inlet - Zero Enrichment (35 min)



## Threonine

- Fragmentation pattern<sup>[7]</sup> defines PSIA targets
- Dual Inlet - Zero Enrichment (24 min)



[6] Shen Y.; et al. J. Dairy Sci., 2014, 97, 679-686.

[7] Zhang P.; et al. Sci. Rep., 2019, 9, 6453.

# Future Frontiers

- PFAS
- Online HPLC methods
- Glucose

Nils Kuhlbusch



# Comprehensive IRMS of sulfates

- Dual Inlet experiment:

Reference: S-3477 (Working Standard)

Sample: S-MIF-1, S-MIF-2


- Reported isotope data:[3]

S-3477 $\delta^{xx}$ [‰]		S-MIF1 $\delta^{xx}$ [‰]		S-MIF2 $\delta^{xx}$ [‰]	
33S	(3.33)	33S	14.81	33S	22.42
34S	6.48	34S	10.26	34S	21.53
36S	(12.35)	36S	19.47	36S	40.73
17O	(6.77)	17O	na	17O	na
18O	13.02	18O	na	18O	na
$\Delta 33$	na (0)	$\Delta 33$	9.54	$\Delta 33$	11.39
$\Delta 36$	na (0)	$\Delta 36$	-0.14	$\Delta 36$	-0.33
$\Delta 17$	na (0)	$\Delta 17$	3.3	$\Delta 17$	3.3

Böhlke, USGS

Samples provided by Joel Savarino  
Institut des Géosciences de  
l'Environnement/CNRS



- Measured Isotope Data:



[3] Geng L.; et al. J. Anal. At. Spectrom., 2019, 34, 1263–1271.

# Comprehensive IRMS of methanesulfonic acid (MSA)

- Easy access to test new hypotheses




[4] Ishino S.; et al. J. Geophys. Res. Atmos., 2021, 126 (6).

[5] Wang K.; et al. Atmos. Chem. Phys., 2021, 21, 8357–8376.

# Analysis of two MSA samples

‘noM0’ experiments:

- **Dual Inlet analysis** of two MSA samples vs a working standard (Sigma-Aldrich)
  - 3x70 min runs → error bars - st. dev. of the 3 runs



Samples provided by Shohei Hattori  
Tokyo Tech & Nanjing

Preliminary results - unpublished data

# Groundwater Sample Table

Table S2. Sample descriptions and isotopic data (non-Orbitrap).

Lab ID	Sample name	Sample type	NO <sub>3</sub> - δ15N per mil	NO <sub>3</sub> - δ18O per mil	NO <sub>3</sub> - δ17O per mil	NO <sub>3</sub> - Δ17O per mil	NO <sub>3</sub> - mM	Cl-/ M/M	SO <sub>4</sub> =/ M/M	Sample description
<b>Reference materials</b>										
USGS32	USGS32	KNO <sub>3</sub> ref.	<b>180.0</b>	25.3	13.0	-0.2	n/a	n/a	n/a	Bohlke et al. (1993)
USGS34	USGS34	KNO <sub>3</sub> ref.	<b>-1.8</b>	<b>-27.8</b>	-14.8	<b>-0.3</b>	n/a	n/a	n/a	Bohlke et al. (2003)
USGS35	USGS35	NaNO <sub>3</sub> ref.	2.7	<b>56.8</b>	51.1	<b>21.6</b>	n/a	n/a	n/a	Bohlke et al. (2003)
N11	N11	KNO <sub>3</sub> ref.	3.6	26.3	13.5	<b>-0.2</b>	n/a	n/a	n/a	
<b>Desert salts</b>										
N-14362	Atacama	soil leachate	0.0	55.1	49.8	21.2	43.9	1.87	1.53	Jackson et al. (2010)
N-13315	Death Valley	soil leachate	3.5	24.2	20.7	8.1	56.6	n/a	n/a	Jackson et al. (2010)
N-14685	Death Valley	soil leachate	-0.8	34.3	30.8	13.0	23.7	13.9	2.29	Jackson et al. (2010)
N-16854	Namibia	brine	9.1	28.6	21.0	6.1	35.8	103	2.05	Jackson et al. (2015)
N-17499	Antarctica	soil leachate	-16.2	82.5	74.5	31.6	86.4	n/a	n/a	Jackson et al. (2015)
N-13825	UAE sabkha	brine	10.6	32.2	25.5	8.7	173.0	40.3	0.02	Jackson et al. (2015)
<b>Fresh groundwater</b>										
N-13993	New Mexico	groundwater	6.8	0.7	0.4	0.1	0.52	33.1	20.1	Jackson et al. (2010)
N-17008	California	groundwater	7.7	3.6	3.3	1.4	1.04	0.36	0.22	Izbicki et al. (2015)
N-13819	New York	groundwater	7.3	3.5	1.8	0.0	0.50	0.63	0.50	Bohlke et al. (2009)
N-10177	Maryland	groundwater	2.7	1.5	0.8	0.0	1.26	0.30	0.19	Bohlke and Denver (1995)
<b>Reference materials as AgNO<sub>3</sub></b>										
N-15038-Ag	USGS34	AgNO <sub>3</sub> ref.	-2.0	-27.9	-14.9	-0.3	n/a	n/a	n/a	Bohlke et al. (2003)
N-15039-Ag	USGS35	AgNO <sub>3</sub> ref.	2.5	56.0	50.7	21.6	n/a	n/a	n/a	Bohlke et al. (2003)
N-15037-Ag	N11	AgNO <sub>3</sub> ref.	3.3	25.6	13.1	-0.2	n/a	n/a	n/a	
<b>Desert salts as AgNO<sub>3</sub></b>										
N-14362-Ag	Atacama	AgNO <sub>3</sub>	-0.2	54.3	49.4	21.2	n/a	n/a	n/a	Jackson et al. (2010)
N-13315-Ag	Death Valley	AgNO <sub>3</sub>	2.9	23.1	20.1	8.1	n/a	n/a	n/a	Jackson et al. (2010)
N-14685-Ag	Death Valley	AgNO <sub>3</sub>	-0.6	33.0	30.2	13.0	n/a	n/a	n/a	Jackson et al. (2010)
N-16854-Ag	Namibia	AgNO <sub>3</sub>	8.8	27.2	20.3	6.1	n/a	n/a	n/a	Jackson et al. (2015)
N-17499-Ag	Antarctica	AgNO <sub>3</sub>	-16.4	80.8	73.6	31.6	n/a	n/a	n/a	Jackson et al. (2015)
N-13825-Ag	UAE sabkha	AgNO <sub>3</sub>	10.3	31.3	25.0	8.7	n/a	n/a	n/a	Jackson et al. (2015)
<b>Fresh groundwater as AgNO<sub>3</sub></b>										
N-13993-Ag	New Mexico	AgNO <sub>3</sub>	6.5	0.7	0.5	0.1	n/a	n/a	n/a	Jackson et al. (2010)
N-17008-Ag	California	AgNO <sub>3</sub>	7.6	3.3	3.1	1.4	n/a	n/a	n/a	Izbicki et al. (2015)

[2] Hilkert, A.; et al. Anal. Chem. 2021, 93, 9139–9148, SI Table.


# 10 reasons to use Orbitrap Exploris Isotope Solutions

## 10 reasons to use Orbitrap Exploris Isotope Solutions


Expand your isotope research

Intramolecular isotopic information is important as it can provide a better understanding of the processes of molecule formation and provide new geochemical proxies for understanding the world around us. Current technology for deriving intramolecular isotopic information is restricted either by large sample sizes, long analysis times or limited applicability across a range of sample types. The powerful, new Thermo Scientific™ Orbitrap Exploris™ Isotope Solutions for natural abundance isotope ratio analysis opens new dimensions for deriving intramolecular isotopic information. Explore why!


### Orbitrap Exploris Isotope Solutions



Thermo Scientific™ Orbitrap Exploris™ 120 MS: the ultimate workhorse for isotope ratio analysis of oxyanions and small organic molecules (<100 m/z)



Thermo Scientific™ Orbitrap Exploris™ 240 MS: the instrument of choice for isotope ratio analysis of oxyanions and medium sized organic molecules (<200 m/z)



Thermo Scientific™ Orbitrap Exploris™ 480 MS: for customers wanting to push the boundaries of isotope research of oxyanions and large organic molecules such as metabolites (<200 m/z)

Learn more at [thermofisher.com/orbitrap-for-isotopes](https://thermofisher.com/orbitrap-for-isotopes)

For Research Use Only. Not for use in diagnostic procedures. © 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. IN001436 0922S

**1** **10x more productivity,** faster analysis by eliminating several sample introduction systems and multiple IRMS systems from the workflow

**2** **10x less sample preparation,** e.g. by avoiding time consuming (and error-inducing) microbial processing, wet chemistry and mineral precipitation

**3** **Reduced sample sizes** due to increased sensitivity compared to classical IRMS approaches

**4** **Direct measurements of individual analytes in liquid samples** without chemical derivatization or conversion

**5** **Analysis of intact molecules** using soft electrospray ionization

**6** **Simultaneous acquisition of all major and some minor isotopologues** by high-resolution accurate mass IRMS

**7** **Extraction of accurate isotopic information from singly- and multi-substituted isotopologues** utilizing methodology for amplifying signals from minor isotopologues

**8** **Position specific isotope information from functional groups** through controlled fragmentation of molecular ions

**9** **Utilization of proven Dual Inlet principles** for determination of accurate isotope ratios of unknown samples relative to a reference

**10** **Automation of sample introduction** by coupling the Orbitrap Exploris Isotope Solutions to Thermo Scientific™ Vanquish™ Neo UHPLC system

thermo scientific